MIRAPONT®

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Instructions for use

Mode d'emploi

MIRAPONT®

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1. Anwendung

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Instructions for use

1. Use

1. Always store Mirapont in a refrigerator.
2. It is recommended to work with gloves.
3. After opening the metal containers with components A and B, stir the respective content strongly and evenly with a stable plastic spatula. The spatulas must not be interchanged. After that, fill the components into the plastic bottles. Shake both components well before use!
4. Excess material at the impression model can be removed with Mirapont Silifit (material saving).
5. Fill the desired quantities of Mirapont A and B into a measuring cup (mixing ratio 1 : 1).

Number of posts	needed quantity component A	component B
4	7 ml	7 ml
8	12 ml	12 ml
12	15 ml	15 ml

6. Fill the measured components A and B into the measuring cup and stir for 15 – 30 seconds until a homogenous compound has been obtained.
7. Dry the impression well before casting. Pour in the Mirapont material in a thin stream from a height of approx. 20 cm. This can also be done with a vibrating device at very low vibration. Too high vibration causes air bubbles.
8. The material can be removed after 2 hours, it can be sawed and ground without problems now.
9. Mirapont has to be isolated against wax with Mirapont Agent Plus.
10. For impressions from polyether or polysulfide (Impregnum, Permadyne, Permplastic, etc.) it is absolutely necessary to isolate the impression surface with the special isolation agent Mirapont Agent Plus prior to casting. It has to be sprayed on from a distance of approx. 20 cm without leaving a dense film. Spray once again directly before casting.
11. Shrinkage of Mirapont is below 0.02 %.
12. Special indications:

Posts for primary telescope- respectively cone crowns, inlays, implants, especially fragile posts, model cast, etc.

When using pin systems, it is advisable to cast the polyurethane prior to inserting the pins into the post in order to avoid air bubbles on the post. Leave the model for 2 hours before removing the impression. Lumps, which might build in the material, do not change its technical characteristics.

2. Characteristics and storage

Mirapont guarantees the precision of your prosthetic work by special fillers in the material. The components A and B are vacuum-packed in metal containers. The scope of supply includes two plastic bottles for preserving the material after opening. Don't store original packing above 20° C for a longer time. Cap the plastic bottles well after each use and keep in a refrigerator. If the products are used regularly, it is sufficient to shake the components A and B well before use, so that the fillers intermingle again. In case the product is not used for 1 – 2 months after opening, shaking is not sufficient. The fillers deposit and have to be intermingled again. Stir the compound with a stable spatula until an even color has been obtained. Note taking a separate spatula for each component, otherwise there are already reactions and the material cures. Mirapont is compatible with the following impression materials: Silicones, alginates, reversible hydrocolloids, polyvinylsiloxanes, thermoplastic materials. It is absolutely necessary to isolate impressions from polyether or polysulfides (Impregnum, Permadyne, Permplastic, etc.) with a special isolation agent (Mirapont Agent Plus) in beforehand. Spray the isolation agent from a distance of approx. 20 cm without leaving a dense film. Dry and spray once more directly before casting. You should wear gloves when working with Mirapont as this is a chemical product, which might affect the skin. After opening the metal containers for the first time, the content should be stirred with an inflexible and stable spatula until a homogenous compound has been obtained. The components can be filled into the plastic bottle as soon as this homogenous compound has been obtained in one color.

It is recommended to use the material within 3 months after opening the metal container. The product should always be stored in a refrigerator.

3. Please note

Shake the bottles again prior to use as the component fillers have deposited at the bottom. The two components should be mixed in equal shares. Ideally you use the provided measuring cups with indicator scale. Mixing time is 15 – 30 seconds. Make sure that casting is done in an even stream (by means of a spatula) or cast on a vibrating unit in order to avoid air bubbles. referably use regular monomer – no acetone – for cleaning of the posts. Cured Mirapont material resists boiling water, steam ejectors and can be washed under hot running water without problems. The material also withstands a temperature of 120° C in the Ivomat device, where it gets a little soft, however, does not deform. For regaining of hardness, cool with an air flow for a short time.

If component B seems to be too hard, the reason may be:

- The bottle has not been capped directly after use so the content became humid
When the material expands and gets porous:
- It got into contact with water before curing

Absolutely avoid contact of the material with water before curing. Cap the bottle directly after use and don't store near steam ejectors. In case the non-hardened material comes into contact with water, the freshly cast impressions should be kept in a pressure container for approx. 15 minutes in order to avoid porosities.

Technical data

Chemical composition:

Component A:	polyol
Component B:	isocyanates
Mixing ratio:	1 : 1
Volume ratio:	50 % A : 50 % B
Mixing time:	20 – 30 seconds
Flow characteristics:	low flow consistency
Processing time:	approx. 3 – 4 minutes
Coagulation:	approx. 10 minutes
Separation from impression:	after 2 hours
Complete curing time:	2 hours, medium – 70 SH D
No chipping and shivering when grinding the material	
Heat resistance:	approx. 120° C
Shrinkage:	0.02 %
(based on a linear casting mold of 10 x 0.5 x 2 cm)	
Contraction not traceable	
Color:	ivory
Storage life of the product	
Unopened:	2 years
Opened:	3 – 4 months (when stored in a cool and dry place)

Mirapont Set incl. Mirapont Agent Plus, Mirapont Silifit, measuring cup, spatula, gloves.

Mirapont 500 g	REF 203 010
Mirapont 1000 g	REF 203 011
Mirapont Agent Plus (200 ml)	REF 203 016
Mirapont Silifit (250 g)	REF 203 015

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Häger & Werken GmbH & Co. KG